



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Men Comfort Knit Red Flame-Resistant Henley

**Product Code :**  
LBNY-00017-11000902



**Description :**

Men Comfort Knit Flame Resistant Henley compliant and a casual, comfortable and stylish part of your FR system.  
One concealed chest pocket with pencil stall and one concealed pencil stall sleeve pocket  
Hidden rib knit cuffs to keep sleeves in place  
With raglan sleeves and underarm gussets for added mobility, so you can bend, reach, stretch and work better  
Fabric: Exclusively with Flame resistant, 5.3 oz. (180 g/m2)  
Blend: 69% Cotton/25% Polyester/6% Polyoxadiazole  
Closure: Four-button placket  
Moisture-wicking, breathable fabric  
Hazard: Flash Fire, Arc Flash  
Wash Care: Home Wash.

**Technical Specification :**

Men Comfort Knit Red Flame-Resistant Henley



**Laboratory Instrument India**

