



Email : sales@laboratoryinstrumentindia.com

Product Name :
Men Comfort Knit Red Flame-Resistant Henley

Product Code :
LBNY-00017-11000902



Description :

Men Comfort Knit Flame Resistant Henley compliant and a casual, comfortable and stylish part of your FR system.
One concealed chest pocket with pencil stall and one concealed pencil stall sleeve pocket
Hidden rib knit cuffs to keep sleeves in place
With raglan sleeves and underarm gussets for added mobility, so you can bend, reach, stretch and work better
Fabric: Exclusively with Flame resistant, 5.3 oz. (180 g/m2)
Blend: 69% Cotton/25% Polyester/6% Polyoxadiazole
Closure: Four-button placket
Moisture-wicking, breathable fabric
Hazard: Flash Fire, Arc Flash
Wash Care: Home Wash.

Technical Specification :

Men Comfort Knit Red Flame-Resistant Henley



Laboratory Instrument India

