



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Men Midweight Comfort Touch Deluxe Insulated Bib Overall

**Product Code :**

LBNY-00017-11000903

**Description :**

Men Midweight Deluxe Insulated Bib Overall offers unbeatable warmth with CAT 2 AR/FR protection and deluxe features like rugged, double-reinforced knees. Hook and loop leg tabs ensure your coverall fits snugly around your boots.

Large patch pockets and bonus utility pockets ensure there's a place for all of your tools.

Features: Reinforced double knees, Elasticized shoulder straps with non-conductive adjustable hardware, Two large center front pockets with pencil stall

For cold jobs, this midweight insulated bib overall offers unbeatable warmth with CAT 3 AR/FR protection and deluxe features like rugged, double-reinforced knees.

Hazard: Flash Fire, Arc Flash

Primary Closure: Concealed Zipper

Closure: Concealed Nomex™ taped brass zipper front

Pocket: Two large front bib pockets with pencil stall, Two large hip pockets, One rule pocket on right leg, Two front pockets

Wash Care: Industrial Laundry - Light Soil, Home Wash

Fabric: Shell: Flame resistant, 7 oz. Twill, Water-repellent, Insulation Package: 14 oz.

Blend: Shell: 88% Cotton/ 12% Nylon, Insulation Package: 50% Modacrylic 50% Cotton

Gender: Male

Materials: Cotton/Nylon

Protection: ATPV Arc Rating: 36

Fit Details: Designed to be worn over other clothing. Fit over work boots.

**Technical Specification :**

---

Men Midweight Comfort Touch Deluxe Insulated Bib Overall



**Laboratory Instrument India**