



Email : sales@laboratoryinstrumentindia.com

Product Name :
Women Royal FR Lab Coats, Knit Cuffs

Product Code :
LBNY-00017-11000907



Description :

Women flame-resistant lab coat with exterior pockets and knit cuffs.
Flame resistant women lab coat is designed for functionality, comfort, and safety
Concealed and reinforced snap front closures for easy donning and doffing
Side vents at hem for ease of walking
Four exterior pockets that provide ample storage
Wash Care : Industrial Laundry - Light Soil
Fabric: Body: 4.5 oz. (150 g/m²) Nomex IIIA Flame Resistant Plain Weave
Blend: Body: 93% Nomex Aramid / 5% Aramid Fibers / 2% Other Fibers
Protection: ATPV Arc Rating: 4.4.
CAT 1 arc rating
Rib knit cuffs
Hazard : Flash Fire, Arc Flash.

Technical Specification :

Women Royal FR Lab Coats, Knit Cuffs



Laboratory Instrument India

