



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Mens Royal Blue Lapel Counter Coats

**Product Code :**

LBNY-00017-11000909



**Description :**

Lab coat with plenty of pocket space between a left chest pocket and two lower ones so you can keep tools handy.

The mens three button front closure on this coat features non-yellowing buttons.

Mens poly cotton blazer-style lapel counter coat with exterior pockets offers a professional appearance.

Fabric: 5 oz. Performance Blend Poplin

Blend: 80% Polyester / 20% Combed Cotton

Primary Closure: Button

Wash Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash

Length: KP10 - 30 in. (based on size medium - back)

Pocket: Left exterior chest pocket, Bartacked pencil stall, Two lower exterior pockets

Finish: Wrinkle Resistant

Closure: Three buttons

Collar: One-piece, lined, notched lapel

One chest pocket and two lower pockets provide storage for tools and personal items

Industrial laundry friendly

Non-yellowing buttons for long-lasting whiteness and professional appearance

Requires minimum ironing thanks to a wrinkle-resistant finish

Blazer styling provides a professional look

Fit Details: Fits true to size. Meant to be worn over clothes.

Gender: Male

Materials: Poplin, Cotton/Polyester.

**Technical Specification :**

---

Mens Royal Blue Lapel Counter Coats



**Laboratory Instrument India**