



Email : sales@laboratoryinstrumentindia.com

Product Name :

FR Straight Fit Jean with Stretch Length: 32 in.

Product Code :

LBNY-00017-11000912



Description :

Flame resistant straight fit jean with length 32 in.

Five jean-style pockets

Felled inseam and seat seam

Soft, pre-washed and sanded fabric

Fabric: Flame-resistant, 12 oz. (405 g/m²) 94% Cotton/ 5% Polyester/ 1% Spandex

Care: Industrial Laundry -Light Soil, Home Wash

One-piece waistband with metal button closure, metal is not exposed to skin

Waistband ticket and size strip

Protection: Arc Rating ATPV 15 calories/cm².

Technical Specification :

FR Straight Fit Jean with Stretch Length: 32 in.



Laboratory Instrument India