



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Men White Gripper-Front Lab Coats

**Product Code :**  
LBNY-00017-11000914



#### **Description :**

Men poly cotton gripper-front lab coat features convenient snaps and exterior pockets.

Primary Closure: Snap/Gripper

Fabric: 5 oz. Performance Blend Cotton

Blend: 80% Polyester / 20% Cotton

Finish: Wrinkle Resistant

Wash Care : Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash

Closure: Four grippers

Collar: One-piece, notched lapel

No buttons mean fewer repairs and additional protection

Requires minimum ironing thanks to a wrinkle-resistant finish

Length: 41.5 in. (based on size medium - back)

Pocket: Left interior chest pocket, two lower exterior pockets

For ultimate protection from leaks and spills in the lab. Made out of Performance Blend poplin and has a buttonless four-gripper closure to keep everything nice and tidy. This lab coat also comes with three large pockets so that you can keep everything you need close by.

Interior chest pocket and two lower pockets provide storage for tools and personal items

Fit Details: Fits true to size. Meant to be worn over clothes.

Side vent openings for easy access to underneath attire

Industrial laundry friendly

Gender: Male.

#### **Technical Specification :**

---

Men White Gripper-Front Lab Coats



**Laboratory Instrument India**