



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

40 Cal Hood Lift FRT 9 oz. Navy 2 PLY Ultrasoft

**Product Code :**

LBNY-00017-11000916



**Description :**

Allows for quick access to fresh air and a clearer view  
Fabric/Material: 9 oz. Ultrasoft Sateen, Double Layer  
Lift front shield makes it easy to obtain fresh air while out of the hazard zone  
Knit gussets at side provide greater range of motion and flexibility  
Premium anti-fog and scratch resistant coating on face shield  
Arc Rating: 40 cal/cm<sup>2</sup>  
CAT 4  
Flash light with mounting clip included  
Vents provide fresh air and increased auditory capabilities  
Pure View face shield technology features a greater range of visibility  
Hard hat included  
Color: Navy  
One size fits most.

**Technical Specification :**

40 Cal Hood Lift FRT 9 oz. Navy 2 PLY Ultrasoft

---

# Laboratory Instrument India

