



Email : sales@laboratoryinstrumentindia.com

Product Name :

AAMI Level 3 Blue Heavyweight Isolation Gowns, Full Back

Product Code :

LBNY-00017-11000917



Description :

AAMI Level 3 isolation gown constructed of a heavy SMS material provides moderate fluid protection.
Strong and breathable
Ideal for extended wear
SMS fabric has inherent low-linting characteristics
No tape tabs or ties to fasten at the neck
Hybrid elastic wrist and thumb hook.
AAMI Level 3
Not made with natural rubber latex
Over-the-head neck style.

Technical Specification :

AAMI Level 3 Blue Heavyweight Isolation Gowns, Full Back



Laboratory Instrument India