



Email : sales@laboratoryinstrumentindia.com

Product Name :
Chemical Resistant FR Lab Coats with Pockets

Product Code :
LBNY-00017-11000919



Description :

Navy chemical resistant lab coat made of flame-resistant CR/FR fabric for multi-hazard protection. Pockets feature covered flaps to secure contents.
Side pass-through allows access to pockets underneath the coat
Durable accessory loop for safety glasses, pens, badges, etc.
Fabric/Material 4.4 oz. CR/FR Aramid Blend FR Nomex knit wrists keep sleeves in place
Bottom side vents for increased mobility
Durable hi-temp plastic snaps for easy donning and doffing
Highly breathable
Easy care and laundering.
Chest pocket with covered flap and pen/pencil slot
Lower pockets with covered flaps to secure contents.

Technical Specification :

Chemical Resistant FR Lab Coats with Pockets



Laboratory Instrument India

