



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Chemical Resistant FR Lab Coats with Pockets

**Product Code :**  
LBNY-00017-11000919



**Description :**

Navy chemical resistant lab coat made of flame-resistant CR/FR fabric for multi-hazard protection. Pockets feature covered flaps to secure contents.  
Side pass-through allows access to pockets underneath the coat  
Durable accessory loop for safety glasses, pens, badges, etc.  
Fabric/Material 4.4 oz. CR/FR Aramid Blend FR Nomex knit wrists keep sleeves in place  
Bottom side vents for increased mobility  
Durable hi-temp plastic snaps for easy donning and doffing  
Highly breathable  
Easy care and laundering.  
Chest pocket with covered flap and pen/pencil slot  
Lower pockets with covered flaps to secure contents.

**Technical Specification :**

Chemical Resistant FR Lab Coats with Pockets



**Laboratory Instrument India**

