



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Chemical Resistant FR Lab Coats, No Pockets

**Product Code :**  
LBNY-00017-11000920



**Description :**

Fabric/Material 4.4 oz. CR/FR Aramid Blend FR Nomex knit wrists keep sleeves in place  
Bottom side vents for increased mobility  
Durable hi-temp plastic snaps for easy donning and doffing  
Chemical Resistant Lab Coat made of flame-resistant CR/FR fabric for multi-hazard protection without weighing you down.  
Highly breathable  
Easy care and laundering.  
Side pass-through allows access to pockets underneath the coat  
Durable accessory loop for safety glasses, pens, badges, etc.

**Technical Specification :**

Chemical Resistant FR Lab Coats, No Pockets



**Laboratory Instrument India**