



Email : sales@laboratoryinstrumentindia.com

Product Name :
Chemical Resistant FR Lab Coats, No Pockets

Product Code :
LBNY-00017-11000920



Description :

Fabric/Material 4.4 oz. CR/FR Aramid Blend FR Nomex knit wrists keep sleeves in place
Bottom side vents for increased mobility
Durable hi-temp plastic snaps for easy donning and doffing
Chemical Resistant Lab Coat made of flame-resistant CR/FR fabric for multi-hazard protection without weighing you down.
Highly breathable
Easy care and laundering.
Side pass-through allows access to pockets underneath the coat
Durable accessory loop for safety glasses, pens, badges, etc.

Technical Specification :

Chemical Resistant FR Lab Coats, No Pockets



Laboratory Instrument India