



Email : sales@laboratoryinstrumentindia.com

Product Name :
Resistance Box

Product Code :
LBNY-0003-1010051



Description :

Resistance Box

Technical Specification :

For accurate work, in moulded box, ohms mentioned on panel.

Split brass contact blocks hold precision cut, interchangeable plug having moulded black fluted tops.

Coils of constantan wire, double silk covered are wound co-efficient.

Range: 500 ?

Accuracy $\pm 0.1\%$.

In bakelite cases, brass blocks, plug type.

Dimensions 200 x 100 x 60mm



Laboratory Instrument India