



Email : sales@laboratoryinstrumentindia.com

Product Name :
Comfort Lightweight Pant Length: 34 in

Product Code :
LBNY-00017-11000924



Description :

Flame resistant lightweight pants
Full crotch gusset
Utility loop extension on front belt loop
Cargo pocket with snap closures on left leg
One-piece waistband with concealed button closure
Two scoop front pockets, concealed pocket on right pocket bag, two patch hip pockets
Fabric: Ripstop Flame-resistant, 6.9 oz. (235 g/m²) 65% Cotton/ 35% Polyester
Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash
Protection: Arc Rating ATPV 9.0 calories/cm²
Double utility pocket on right leg
Double articulated knees one in. silver trim sewn just below reinforced knee.

Technical Specification :

Comfort Lightweight Pant Length: 34 in



Laboratory Instrument India

