



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
200 Totally Encapsulating Suits

**Product Code :**  
LBNY-00017-11000925



**Description :**

Expanded back, rear entry, side air inlet, double storm flap with hook and loop closure, cloth zipper  
For Use With (Application): Ideally used with wider range of chemicals.  
Closure Type: Double Storm Flap with Hook and Loop Seal  
Seam Type: Sealed, Taped  
Resistance: Chemicals  
Back Style: Expanded  
Size: 2X-Large  
Color: Gray

**Technical Specification :**

200 Totally Encapsulating Suits



**Laboratory Instrument India**