



Email : sales@laboratoryinstrumentindia.com

Product Name :
200 Totally Encapsulating Suits

Product Code :
LBNY-00017-11000925



Description :

Expanded back, rear entry, side air inlet, double storm flap with hook and loop closure, cloth zipper
For Use With (Application): Ideally used with wider range of chemicals.
Closure Type: Double Storm Flap with Hook and Loop Seal
Seam Type: Sealed, Taped
Resistance: Chemicals
Back Style: Expanded
Size: 2X-Large
Color: Gray

Technical Specification :

200 Totally Encapsulating Suits



Laboratory Instrument India