



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Hot Rods Winter Liners: Regular Length

**Product Code :**  
LBNY-00017-11000929



**Description :**

Insulated ear barriers provide extra protection against wind, water and cold  
Strip of reflective film above the back of neck offers visibility in low-light or dark hours of the day.  
Feature two built-in pockets to hold regular-size Hot Rods warming packs  
Short elastic strip in the forehead area keeps liner snug to face.

Two layers

Color: Navy

Three layers

100% quilted nylon shell

2mm foam center layer

100% cotton twill top

Fleece lining

Flame-retardant until washed

Fleece lining

Color: Green.

**Technical Specification :**

Hot Rods Winter Liners: Regular Length

---

# Laboratory Instrument India

