



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hot Rods Winter Liners: Regular Length

Product Code :
LBNY-00017-11000929



Description :

Insulated ear barriers provide extra protection against wind, water and cold
Strip of reflective film above the back of neck offers visibility in low-light or dark hours of the day.
Feature two built-in pockets to hold regular-size Hot Rods warming packs
Short elastic strip in the forehead area keeps liner snug to face.

Two layers

Color: Navy

Three layers

100% quilted nylon shell

2mm foam center layer

100% cotton twill top

Fleece lining

Flame-retardant until washed

Fleece lining

Color: Green.

Technical Specification :

Hot Rods Winter Liners: Regular Length

Laboratory Instrument India

