



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Sleeves with Aramid Fiber

**Product Code :**  
LBNY-00017-11000933



**Description :**

Constructed of 100% aramid fiber  
Tubular knit design with optional thumb hole keeps sleeve in place  
Double-ply construction  
Sleeves protect against cuts and abrasions from glass, metals and other sharp objects.  
Cut level 3  
One size fits all.

**Technical Specification :**

Sleeves with Aramid Fiber



**Laboratory Instrument India**