



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bulk Dispensers

Product Code :
LBNY-00017-11000934



Description :

Freestanding units hold everything from bouffant caps to shoe covers— each in its own compartment
Sturdy, see-through acrylic for unobstructed viewing
Open top offers easy access
Flame-treated edges
Rubber pads prevent slipping and scratching
Three-compartment dispensers measure 13L x 18.5W x 7 in. H (33 x 47 x 17.8cm).

Technical Specification :

Bulk Dispensers



Laboratory Instrument India