



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Standard Weight Polyethylene Shoe Cover

**Product Code :**  
LBNY-00017-11000938



**Description :**

Manufactured with 100% cross-linked polyethylene with heat-sealed seams. Strong, durable, impervious to non-hazardous liquids, and critical environment compatible.

Provides full liquid protection and is commonly used in microelectronics and pharmaceutical industries.

This product is recommended for use in a cleanroom Class 100 -1,000 critical environment

All shoe covers are tested and manufactured in under inspection and stringent process control to ensure quality standard and product specification compliance.

**Technical Specification :**

Standard Weight Polyethylene Shoe Cover



**Laboratory Instrument India**