



Email : sales@laboratoryinstrumentindia.com

Product Name :
Standard Weight Polyethylene Shoe Cover

Product Code :
LBNY-00017-11000938



Description :

Manufactured with 100% cross-linked polyethylene with heat-sealed seams. Strong, durable, impervious to non-hazardous liquids, and critical environment compatible.
Provides full liquid protection and is commonly used in microelectronics and pharmaceutical industries.
This product is recommended for use in a cleanroom Class 100 -1,000 critical environment
All shoe covers are tested and manufactured in under inspection and stringent process control to ensure quality standard and product specification compliance.

Technical Specification :

Standard Weight Polyethylene Shoe Cover



Laboratory Instrument India