



Email : sales@laboratoryinstrumentindia.com

Product Name :
CTI Polypropylene Shoe Covers

Product Code :
LBNY-00017-11000939



Description :

Strong and durable shoe covers are heavily textured for a sure grip and are compatible with class 100 environments.
Individually inspected for stitching, sizing, materials and workmanship
Lot traceability
Lightweight material offers comfort and breathability, while maintain snug, full coverage
Color: Blue.

Technical Specification :

CTI Polypropylene Shoe Covers



Laboratory Instrument India