




Email : sales@laboratoryinstrumentindia.com

Product Name :
Locking N95 Mask Dispenser

Product Code :
LBNY-00017-11000944



Description :

Made with extremely durable, break resistant 3/16 in. thick PETG
Choose from keyhole slots for wall mounting, or powerful magnets.
Holds and dispenses one box of N95 Masks.
Hardware not included for wall mount unit (-W). Customer must provide their own screw/anchor for their wall type.
Locking, sliding lid keeps dust out.

Technical Specification :

Locking N95 Mask Dispenser



Laboratory Instrument India