




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Locking N95 Mask Dispenser

**Product Code :**  
LBNY-00017-11000944



**Description :**

Made with extremely durable, break resistant 3/16 in. thick PETG  
Choose from keyhole slots for wall mounting, or powerful magnets.  
Holds and dispenses one box of N95 Masks.  
Hardware not included for wall mount unit (-W). Customer must provide their own screw/anchor for their wall type.  
Locking, sliding lid keeps dust out.

**Technical Specification :**

Locking N95 Mask Dispenser



**Laboratory Instrument India**