




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Shoe Cover Dispenser Type 1

**Product Code :**  
LBNY-00017-11000947



**Description :**

Holds and dispenses a box of shoe covers.  
Durable PETG construction.  
Super durable and cleanable (alcohol resistant)  
Mounts with opening on left, top, or right.  
Choose from mounting methods:  
Powerful Magnets  
Keyhole slots.

**Technical Specification :**

Shoe Cover Dispenser Type 1



**Laboratory Instrument India**