



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Large PPE Box Disp Mount

**Product Code :**  
LBNY-00017-11000948



**Description :**

Holds and dispenses large box of PPE, such as bouffants, shoe covers, or hair nets. Durable, cleanable PETG construction.

Holds and dispenses large box of PPE such as bouffants, shoe covers, hair nets, etc. Choose from magnet mount, VHB tape mount, or wall mount.

Durable, cleanable PETG construction

Choose from VHB tape mount (very high bond tape included), or magnet mount (powerful magnets included).

Magnet mount (powerful magnets included).

VHB tape mount (very high bond tape included)

Wall mount (screws and anchors included).

**Technical Specification :**

Large PPE Box Disp Mount



**Laboratory Instrument India**