



Email : sales@laboratoryinstrumentindia.com

Product Name :
Earloop Mask Dispenser

Product Code :
LBNY-00017-11000949



Description :

The Ear loop Mask Dispenser holds a stack of bulk earloop masks.
Durable, cleanable PETG construction.
Unique design helps to limit mask dispense to one piece each time.
Masks dispense from bottom front.
See image for dimensions.
Choose your mounting method: Wall mount (Screws and anchors included), Magnet Mount (Powerful magnets included), VHB Tape Mount (VHB Tape included), or steel stand mount.

Technical Specification :

Earloop Mask Dispenser



Laboratory Instrument India