




Email : sales@laboratoryinstrumentindia.com

Product Name :
Mask Box Holder with Lock and Key

Product Code :
LBNY-00017-11000950



Description :

Holds a box of ear loop masks, while keeping the contents secure with a key lock.
For stand mount version, choose from sleeve sign, or no sign.
Choose the amount of gloves the glove holder holds.
Choose your mounting method: Stand (steel stand included), Wall mount (screws and anchors included), Magnet Mount (powerful magnets included), or Tape Mount (powerful VHB tape included).
Available holders for 2, 3, or 4 glove capacity.

Technical Specification :

Mask Box Holder with Lock and Key



Laboratory Instrument India