



Email : sales@laboratoryinstrumentindia.com

Product Name :
Protective Wear Zipper Front Reusable Cleanroom Frock

Product Code :
LBNY-00017-11000951



Description :

Reusable frocks with coil zipper enclosure are suitable for a variety of cleanroom applications including pharmaceutical, biotechnology, and medical device manufacturing.

Wide-pull zipper with a single snap at bottom

Cuff Sleeves

Fluid and bacteria resistant Integrity high performance carbon striped fabric

Mandarin collar with two snaps at top

Silicon-free

Laser-cut, pre-serged seams

ESD carbon stripes provide electrostatic dissipation.

Technical Specification :

Protective Wear Zipper Front Reusable Cleanroom Frock



Laboratory Instrument India