



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Protective Wear Taffeta General Purpose Cleanroom Hood

**Product Code :**  
LBNY-00017-11000952



**Description :**

Taffeta hood is ideal for general cleanroom applications.  
Snaps to adjust fit.  
Open Face Style  
Supports detachable face mask.

**Technical Specification :**

Protective Wear Taffeta General Purpose Cleanroom Hood



**Laboratory Instrument India**