



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mens Specialized Lapel Counter Coats

Product Code :
LBNY-00017-11000953



Description :

Professional style lab coat features a square-hem, three gripper front closure, and a lined collar.
Fabric: 5 oz. Performance Blend Poplin
Blend: 80% Polyester/20% Combed Cotton
Finish: Wrinkle Resistant
Primary Closure: Snap/Gripper
Collar: One-piece, lined, notched lapel
Length: KP16 - 30 in. (based on size medium - back), KP17 - 28 in. (based on size medium - back)
Pocket: Two lower exterior pockets
Wash Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash
Closure: Three grippers

Features:

Materials: Poplin, Cotton/Polyester.
Two lower exterior pockets for additional storage
Squared-hem for a professional look
No buttons mean fewer repairs and additional protection
Requires minimum ironing thanks to a wrinkle-resistant finish
Gender: Male
30 in. back length on size medium
Fit Details: Fits true to size. Meant to be worn over clothes.

Technical Specification :

Mens Specialized Lapel Counter Coats



Laboratory Instrument India