



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Mens Specialized Lapel Counter Coats

**Product Code :**  
LBNY-00017-11000953



**Description :**

Professional style lab coat features a square-hem, three gripper front closure, and a lined collar.  
Fabric: 5 oz. Performance Blend Poplin  
Blend: 80% Polyester/20% Combed Cotton  
Finish: Wrinkle Resistant  
Primary Closure: Snap/Gripper  
Collar: One-piece, lined, notched lapel  
Length: KP16 - 30 in. (based on size medium - back), KP17 - 28 in. (based on size medium - back)  
Pocket: Two lower exterior pockets  
Wash Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash  
Closure: Three grippers

**Features:**

Materials: Poplin, Cotton/Polyester.  
Two lower exterior pockets for additional storage  
Squared-hem for a professional look  
No buttons mean fewer repairs and additional protection  
Requires minimum ironing thanks to a wrinkle-resistant finish  
Gender: Male  
30 in. back length on size medium  
Fit Details: Fits true to size. Meant to be worn over clothes.

**Technical Specification :**

---

Mens Specialized Lapel Counter Coats



**Laboratory Instrument India**