



Email : sales@laboratoryinstrumentindia.com

Product Name :

White Biohazard Suit with Elastic Wrists and Ankles

Product Code :

LBNY-00017-11000957



Description :

These disposable coveralls are designed with fabric and seams that have been rigorously tested for protection against blood and blood-borne pathogens.
Zipper front with a sealable storm flap helps keep contaminants from entering or seeping through exposed zipper seams
Red fabric pull tab on the zipper makes it quick and easy to find so you can suit up or down in no time
Durable fabric is built to withstand tears and abrasions in tough working environments
Double-sided tape around wrists helps keeps gloves in place during donning
Storm flap over zipper helps prevent contaminants from entering through otherwise exposed zipper seams
Taped seams help keep contaminants out for added protection, while elastic wrists and ankles keep coveralls secure and in place
Convenient thumb and finger loops help keep sleeves from riding up while you work
Possible applications include emergency response, healthcare, category 3 black water, crime scene remediation, and animal research.
Roomy one-piece raglan sleeves resist tears around shoulders
A gusseted crotch gives the user more freedom to squat, bend, and crouch without ripping fabric.

Technical Specification :

White Biohazard Suit with Elastic Wrists and Ankles

Laboratory Instrument India

