



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Micro Guard MP Microporous Sleeves with Elastic Closures

**Product Code :**

LBNY-00017-11000958



**Description :**

Using state-of-the-art lamination technology, Micro Guard MP microporous 18 in. sleeves combine a soft, non-woven fabric with a microporous film.

Breathable fabric allows moisture to escape, so you're cooler and more comfortable even as you sweat

Lightweight, tear-resistant liquid and particulate barrier helps safeguard from blood penetration, various industrial chemicals, along with dust, dirt, and grime

Secure elastic closures on both ends help keep sleeves in place even with repetitive movements, allowing you to stay on task without having to continually readjust your sleeves throughout your work day

Environmentally friendly; can be recycled after use

Passes testing for blood penetration.

Generous 18 in. length helps cover more of your arms to decrease the risk of exposing skin or underlying clothing to workplace contaminants.

**Technical Specification :**

Micro Guard MP Microporous Sleeves with Elastic Closures



**Laboratory Instrument India**

