



Email : sales@laboratoryinstrumentindia.com

Product Name :
Evaporative Microfiber Cooling Towel

Product Code :
LBNY-00017-11000960

COOLING RELIEF

AS WATER EVAPORATES,
TOWEL CREATES A
COOLING EFFECT THAT
LASTS FOR HOURS



Description :

Microfiber Cooling Towel is designed to create a long-lasting cool feeling while absorbing and holding water. Made with lightweight, absorbent microfiber material
Stays soft either wet or dry
Machine washable with mild detergent, air dry only
Reduces the effects of heat exhaustion and fights fatigue.
Microfiber material hold 30% more water and cools 30% longer
Run under water for one minute or until saturated, wring out excess water for drip-free cooling sensation as it evaporates
UPF 50 provides sun protection from harmful UV rays.

Technical Specification :

Evaporative Microfiber Cooling Towel



Laboratory Instrument India