



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Retort Stand With Clamp

**Product Code :**  
LBNY-0003-1010055



**Description :**

Retort Stand With Clamp

**Technical Specification :**

Steel base with tapped hole 10 x 1.5mm metric thread in centre size 160 x 100mm, metal rod to fit the hole of length 760mm with clamp.  
Finished in blue and fitted with rubber feet.



**Laboratory Instrument India**