



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Type R Class 3 Two-Tone Vest

**Product Code :**  
LBNY-00017-11000961



**Description :**

Mesh fabric, zipper closure  
Material: Polyester  
Color: Lime  
No. of Pockets: 1 inside.  
Size: 4X-Large to 5X-Large.

**Technical Specification :**

Type R Class 3 Two-Tone Vest



**Laboratory Instrument India**