



Email : sales@laboratoryinstrumentindia.com

Product Name :
Class E Two-Tone Pants

Product Code :
LBNY-00017-11000962



Description :

Class E Two-Tone Pants are constructed of cool mesh with wide leg openings and hook and loop closures. Wide leg openings with hook and loop closures
One right leg cargo pocket.
Reflective set on 4.5 in contrasting background
Drawstring waist
Polyester mesh (3.3 oz.)
2 in. ANSI-compliant reflective.

Technical Specification :

Class E Two-Tone Pants



Laboratory Instrument India