



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Type P Class 2 Expandable Public Safety Vest

**Product Code :**  
LBNY-00017-11000964



**Description :**

Vest with adjustable sides.  
Ideal for wearing over turnout gear.  
Material: Polyester.  
Color: Lime  
Size: 3X-Large Plus.

**Technical Specification :**

Type P Class 2 Expandable Public Safety Vest



**Laboratory Instrument India**