



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Type R Class 3 FR Modacrylic Vest

**Product Code :**  
LBNY-00017-11000967



**Description :**

Class 3 Flame Resistant Modacrylic Vest constructed of modacrylic mesh fabric with a hook and loop closure and D-ring access.

Flame Resistant hook and loop closure with D-ring access slot on the back of the vest

Modacrylic mesh 5.6 oz. body and binding (87% Modacrylic and 13% Aramid Fibers)

One interior pocket and one dual pen pocket for additional storage

Single radio mic tab made from FR modacrylic material.

Fabric passes; ATPV of 4.6 cal/cm<sup>2</sup>

2 in. FR Reflective.

**Technical Specification :**

Type R Class 3 FR Modacrylic Vest



**Laboratory Instrument India**