



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Type R Class 2 Two-Tone X-Back Vest

**Product Code :**  
LBNY-00017-11000970



**Description :**

Design with D-ring access and multiple pockets.

Material: Polyester

No. of Pockets: 2 inside, 4 outside.

Color: Orange

Size: 4X-Large to 5X-Large

**Technical Specification :**

Type R Class 2 Two-Tone X-Back Vest



**Laboratory Instrument India**