



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Rheostat

**Product Code :**  
LBNY-0003-1010056



**Description :**

Rheostat

**Technical Specification :**

The tube is carrying upon Bakelite rider and the heavy duty sliding contact is of the phosphor bronze with brass rod 4 mm socket terminals are fitted, allowing for rheostat or potentiometer connection  
Total resistance and maximum current capacity is marked on the slider  
Total Permanent current loading not exceed 75% of the marked values  
Resistance: 50 ohms.  
Current: 2 Amps.  
Resistance sliding contact tube size 200 x 38mm  
Overall size: 255 x 52 x 92mm.



**Laboratory Instrument India**