



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Double Bulk Cone Mask Dispenser

**Product Code :**  
LBNY-00017-11000974



**Description :**

Holds two stacks and/or two boxes of cone masks. Sliding lid allows access to either side.

Material: PETG

Mounting: Magnetic Mount

Dimensions (D x W x H) Interior: 5.5 x 5.5 x 7 in.

Width (English): 11.5 in.

Depth (English): 5.5 in.

**Technical Specification :**

Double Bulk Cone Mask Dispenser



**Laboratory Instrument India**