




Email : sales@laboratoryinstrumentindia.com

Product Name :
Gravity Close Bin for Bulk Items

Product Code :
LBNY-00017-11000975



Description :

Constructed from super durable 1/4 in. thick PETG
Closing door helps to prevent dust accumulation
Holds and dispenses bulk PPE, tubes, supplies, etc.
Auto close bin for bulk items helps to prevent dust accumulation.
Gravity close bin
Align multiple units for different PPE types, sizes, etc.
Choose from magnet mount, VHB Tape mount, or wall mount.

Technical Specification :

Gravity Close Bin for Bulk Items



Laboratory Instrument India