



Email : sales@laboratoryinstrumentindia.com

Product Name :
Polypropylene Lab Coats, Knit Wrists and Collar

Product Code :
LBNY-00017-11000979



Description :

Lab coats with a heavier polypropylene alternative for splash protection from non-hazardous liquids.
Provides non-hazardous splash and dry particulate protection
Snap front
Three total pockets, one breast and two hip
White, knit cuffs, knit collar.
Lab coats made with heavier 1.75 oz polypropylene
Packaged 25 per case.

Technical Specification :

Polypropylene Lab Coats, Knit Wrists and Collar



Laboratory Instrument India