



Email : sales@laboratoryinstrumentindia.com

Product Name :
Polypropylene Lab Coats, Elastic Wrists

Product Code :
LBNY-00017-11000980



Description :

Lab coat offers lightweight splash protection from non-hazardous liquids and dry particulates
Non-woven material lab coats offer lightweight splash protection from non-hazardous liquids and dry particulates.

Snap front

No pockets, elastic wrists.

Packaged 25 per case

1.25 oz spun bond polypropylene.

Technical Specification :

Polypropylene Lab Coats, Elastic Wrists



Laboratory Instrument India