



Email : sales@laboratoryinstrumentindia.com

Product Name :
Specialized Lab Coats

Product Code :
LBNY-00017-11000981



Description :

Wash Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash
Pocketless lab coat features five snaps for secure closure and side vent openings to allow for mobility and comfort.

Features: No buttons mean fewer repairs and additional protection

Requires minimum ironing thanks to a wrinkle-resistant finish

Side vent openings for easy access to underneath attire.

Fabric: 5 oz. Performance Blend Poplin

Primary Closure: Snap/Gripper

Blend: 80% Polyester/20% Combed Cotton

Finish: Wrinkle Resistant

41.5 in. back length on all sizes

Closure: Five grippers

Collar: One-piece, lined, notched lapel

Length: 41.5 in. (all sizes - back)

Pocket: Pocketless

Fit Details: Fits true to size. Meant to be worn over clothes.

Materials: Poplin, Cotton/Polyester.

Gender: Male.

Technical Specification :

Specialized Lab Coats



Laboratory Instrument India