



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Mens Lapel Counter Coats

**Product Code :**  
LBNY-00017-11000982



**Description :**

Mens professional counter coat features spacious pockets for storage of tools and personal items. Blazer style available in sizes x-small through 4x-large.

Wash Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash

Primary Closure: Button

Materials: Poplin, Cotton/Polyester.

Finish: Wrinkle Resistant

Closure: Three buttons

Collar: One-piece, lined, notched lapel

Fabric: 5 oz. Performance Blend Poplin

Blend: 80% Polyester/20% Combed Cotton

Length: 30 in. (based on size medium - back)

28 in. (based on size medium - back)

Features: Non-yellowing buttons for long-lasting whiteness and professional appearance

One chest pocket and two lower pockets provide storage for tools.

30 in. back length on size medium

Fit Details: Fits true to size. Meant to be worn over clothes.

Pocket: Left exterior chest pocket, Barracked pencil stall, Two lower exterior pockets

Requires minimum ironing thanks to a wrinkle-resistant finish

Blazer styling provides a professional look

Gender: Male.

**Technical Specification :**



## Laboratory Instrument India